

Development of a preclinical porcine model to evaluate hernia mesh migration, folding and dislocation


Arthur Jourdan, Anthony Végleur, Amandine Radlovic, Ludovic Boure, Anicet Le Ruyet
 Medtronic, Core surgical Innovations, Trévoux, France

Context

New mesh development

How can we ensure that new implant designs do not promote failure mechanisms?

- Development of a porcine model for the evaluation of mesh performances
- Combination of *in vivo* and *ex vivo* testing - Investigation of 3 failure mechanisms:



Comparison of 3 non-resorbable mesh technologies:

- Lightweight flatsheet (LWF) - Optilene® Mesh LP - 39 g/m²
- Heavyweight flatsheet (HWF) - Bard® mesh - 99 g/m²
- Heavyweight self-fixating mesh (SFM) - ProGrip™ laparoscopic self-fixating mesh - 147 g/m²

Phase 1: *in vivo* testing

Migration and folding and bunching evaluation, one-week post-surgery

M&M

- Partially radio-opaque 13x9 cm mesh with micro-beads coating
- Laparoscopic bilateral mesh placement + ø3cm defect creation
- AW CT-scan at t0 and 1-week + Post-processing

Abdominal wall Viscera Peritoneum

Meshes

Coring device

Meshes


t0 Beads

Height (mm)


Anterior (mm)

RESULTS - one-week post-surgery - n=4 / mesh type

- High migration and folding and bunching of the LWF
- Less migration and folding of the SFM and HWF



- Mesh dislocation: 100% LWF, 75% HWF and 25% SFM
- Seroma: 75% LWF, 100% HWF and 0% SFM



➤ Reduced risk of migration and bunching with SFM

Phase 2: *ex vivo* testing

Dislocation dynamic bench testing

M&M

- Explanted abdominal wall samples
- 250 mmHg cyclic pressure impacts until mesh dislocation^{1,2,3}
- Simulating a critical postoperative situation²

Compressed air inlet Pressure

Pressure chamber

Pressure sensor Plexiglas

13x9 cm Meshes AW sample

Top view

Bottom view


15 cm ø3cm

Cyclic pressure test (Bottom view)

Initial state Cycle 1 Cycle 150 Cycle 400 Dislocation

Pressure cycles

RESULTS



- Ranking was aligned with *in vivo* results
- SFM provided better protection of the repair

- Mechanisms of failure**
 - Mesh creeping³
 - Defect opening
 - Insufficient overlap
 - Mesh dislocation
- 1 case of LWF mesh rupture

